

## BA2000i Microscope

Items	Spec and Parameter	BA2203i	BA2303i	Order No.
Eyepiece	10X-20	● ●	● ●	EW10
ISO Plan Achromatic Objectives	Plan 10X/0.25	●	●	B-OPi10a
	Plan 20X/0.40	●	●	B-OPi20a
	Plan 40X/0.66	●	●	B-OPi40a
	Plan 100X/1.25(Oil)	●	●	B-OPoi100a
Main Body	Quintuple Nosepiece	●	●	A08a
Seidentopf Binocular Head		●		TBR11
Seidentopf Trinocular Head			●	TTR04
Mechanical Stage	Moving Range: 76×50 Size: 160×135	●	●	S26
Condenser Focusing Team	Rack and Pinion Focusing Mechanism	●	●	CA03a
Condenser	N.A. 1.25			C23
Field Diaphragm	Ø 2-ø35	●	●	DI06
Immersion Oil		●	●	
Filter	Ø 32 (Blue)	●	●	CU2-1

## BA2000i Microscope

### BA2000i Series Optional Accessories

Items	Spec and Parameter	Order No.	
Eyepieces	10X-18 (Pointer)	Erg10	
	16X-10	EP16a	
	16X-14	EW16	
IOS Plan Achromatic Objective	Plan 4X/0.10	B-OPi4a	
Phase Contrast Attachment	IOS Plan Positive Ph. Objectives	Plan ph10x/0.25	B-OPip10a
		Plan ph20x/0.40	B-OPip20a
		Plan ph40x/0.66	B-OPip40a
		Plan ph100x/1.25(Oil)	B-OPoip100a
	IOS Plan Negative Ph. Objectives	Plan ph10x/0.25	B-OPin10
		Plan ph20x/0.40	B-OPin20
		Plan ph40x/0.66	B-OPin40
		Plan ph100x/1.25(Oil)	B-Opoin100
	Phase Contrast Condenser		C22
	Centering Eyepiece		ECi
Dark Field Attachment	Dark Field Condenser (N.A. 0.83-0.91)	C07	
Polarizing Attachment	Polarizer Slide	PO01	
	Analyzer Slide	PA01	
Photo Attachment	10X Reticle Eyepiece	EPr10b	
	2.5X Photo Eyepiece	EPHn2.5	
	4X Photo Eyepiece	EPHn4	
	6.3X Photo Eyepiece	EPHn6.3	
	135 Photo Device (With 4x Photo Eyepiece)	MP03	
TV and Video Attachment	C-mount Adaptor、Color or Black and White Monitor、CCD Camera	TV02	
Epi-illumination Attachment	6V20W Halogen Tungsten Lamp (With Reflector)	H09	
Side Viewing Attachment	Two View Field Space 600mm	TB06	
Projection Attachment		PJ04	
Filter	Ø 32 (Green、Amber)	CU02-2、3	
Image Analysis and Disposal System			

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139