

XSZ-4Ga microscope

Model		4Ga101	4Ga201	4Ga301	Order No.
Items Spec. and parameter					
Eyepiece	Eyepiece WF10X/.18	●	●●	●●	EP10i
	Eyepiece WF16X /.10	○	○○	○○	EP16b
Objective	Achromatic Objective 4X/0.10	●	●	●	B-OA4d
	Achromatic Objective 10X/0.25	●	●	●	B-OA10e
	Achromatic Objective 40X/0.65	●	●	●	B-OSP40f
	Achromatic Objective 100X/1.25 (oil)	●	●	●	B-OAo100d
Head	Monocular head (Inclined 30° , Rotatable 360°)	●			62X.1
	Seidentopf Binocular head (Inclined 30° , Rotatable 360° , Interpupillary Distance: 55-75mm)		●		62X.2
	Seidentopf Trinocular head (Inclined 30° , Rotatable 360° , Interpupillary Distance: 55-75mm)			●	62X.7
	Teaching head (Inclined 30° , Rotatable 360°)	○			62X.8
Main body		●	●	●	62X.3
Stage	Large scale stage: moving range: 78×30mm	○	●	●	
	Ordinary stage: moving range: 78×30mm	●			
Abbe condenser:		●	●	●	
Base	Ordinary base	●	●	●	62X.5
	Chargeable base (for areas usually lack of electricity)	○	○	○	62X.5A

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

CE ISO9001
ISO14001



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139