

I-L3 Serials Microscope

Model	I-L3
Eyepiece	Wide field WF10X(ϕ 22mm)
Objective	Infinity plan achromatic objectives with long working distance (no cover glass)
	PL L5X/0.12(Work distance): 26.1 mm
	PL L10X/0.25(Work distance): 20.2 mm
	PL L40X/0.60 ∞ Spring ∞ (Work distance): 3.98 mm
	PL L60X/0.70 ∞ Spring ∞ (Work distance): 3.18 mm
Eyepieces tube	Trinocular inclined 30°, can be shot in 100% light flux.
Epi-illumination system	6V 30W, halogen lamp, adjustable brightness
	Integrated field diaphragm, aperture diaphragm and filter switching device. Push-pull type analyzer and polarizer
Focus system	Coaxial coarse/fine focus system, with tension adjustable and limit stopper, minimum division of fine focusing: 2 μ m.
Nosepiece	Quintuple (Backward ball bearing inner locating)
Stage	Double layer mechanical (Size:210mmX140mm,movingrange:75mmX50mm)
Transmitted illumination system	Abbe condenser NA.1.25 Rack & pinion adjustable
	Blue filter and Ground glass
	Collector for illumination with halogen lamp (With field diaphragm)
	6V 30W halogen lamp, adjustable brightness



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

I-L3 Serials Microscope

Optional accessories

Name	Sort/Technique parameter	NO.
Eyepiece	Dividing eyepiece(field number: μ 22mm) 0.10mm/Div	1121012
Objective	Infinity plan achromatic objectives with long working distance (no cover glass)	
	PL L 20X/0.40(Work distance): 8.80mm	2830120
	PL L50X/0.70(Work distance): 3.68 mm	2830150
	PL L80X/0.80(Work distance): 1.25 mm	2830180
	PL L 100X/0.85(Work distance):0.4mm	2250111
Filter	Green filter	115002
	Yellow filter	115003
CCD adapter	0.4X	810001
	0.5X	810004
	1X	810002
	0.5X with dividing 0.1mm/Div	810003
Camera	DV-1 Video output(380/520 TV line) USB output 0.42 M pixel	800001
	DV-2 With USB output 1.3M,2.0M,3.0M pixel	800003
	DV-3 With video output(380/520 TV line)	800005
Digital camera adapter	CANON(A570,A610,A620,A630,A640,A650,EF) NIKON(F)	820001



Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139