

XDS-1 Inverted Biological Microscope

Order No.	Items	Standard Outfits
	Main Body (including left & right supporting arm plate, 4 piece of Allen screws & a wrench)	●
EW10	Eyepiece WF10×/ø20	● ●
EW16	Eyepiece WF16×/ø14	● ●
EC11	Centering Eyepiece DZ11×	●
EPr10b	Plan Article Eyepiece PF10×	●
B-OLPr10	LWD Plan Achromatic Objective 10/0.25,160/1.5,W.D.=7.9mm	●
B-OLPr25	LWD Plan Achromatic Objective 25/0.40,160/1.5,W.D.=5mm	●
B-OLPr40	LWD Plan Achromatic Objective 40/0.60,160/1.5,W.D.=3mm	●
B-OLPnr10	LWD Plan Phase Contrast Objective 10/0.25,160/1.5,W.D.=7.9mm	○
B-OLPnr25	LWD Plan Phase Contrast Objective 10/0.40,160/1.5,W.D.=5mm	●
B-OLPnr40	LWD Plan Phase Contrast Objective 40/0.60,160/1.5,W.D.=3mm	○
	35mm DF Camera With Shutter Control Line	○
33X-4	Set for Stage Insert Plate	●
9J-8	Stage Insert plate (III)	●
9J-9	Stage Insert plate (IV)	●
33X-5	Stage Insert plate (V)	●
CL02	Stage Clips	● ●
7X.5.1	Cultural Flask (bottom thickness 1.5mm)	● ●
7X.5.2	Cultural Vessel (bottom thickness 1.5mm)	● ●
RP06	Annular Spot 10×	○
12X.2.3-19,20	Annular Spot 25×	●
RP07	Annular Spot 40×	○
	Halogen Tungsten Lamp 6V20W (G4)	● ● ● ●
9J.7	Photo Adaptor With Photo Eyepiece S5×	●
33X.3.2.1	Movable Condenser Insert Plate	●
TBR03	Seidentopf Binocular Head (inclined 45°)	●
	Micrometer (division 0.01mm)	○
EPHp6.3	Photo Eyepiece S6.3×	○
9J.15	Photo Eyepiece S0.65×For CCD	○
CU43-3	C-mount	○
	CCD TV & Record Attachment	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139