

XJP-6A Serials Microscope

Order No.	Item	Standard Outfits	
		XJP-6	XJP-6A
	Main Body	●	●
S09	Mechanical stage with a stroke of 75 by 50 and size of 180 by 155	●	●
EP10	Plan eyepiece 10X-18	●●	●●
EP12.5	Plan eyepiece 12.X-14	●●	●●
EH5	Huygenian ocular 5X-20	○○	○○
EW10	Plan eyepiece 10X-20	○○	○○
Epr10b	Plan reticle eyepiece 10X-14	●	●
B-OAr4	Achromatic objective 4/0.10	●	
M-OAr10	Achromatic objective 10/0.25	●	
M-OSPr40	S-plan achromatic objective 40/0.65	●	
M-OAor100	Achromatic objective 100/1.25 oil	●	
B-OPr2.5	Plan achromatic objective 2.5/0.07	○	○
B-OPr4	Plan achromatic objective 4/0.10		○
M-OPr10	Plan achromatic objective 10/0.25		●
M-OPr20	Plan achromatic objective 20/0.35	○	●
M-OPr40a	Plan achromatic objective 40/0.65		●
M-OSPor100	S-plan achromatic objective 100/1.25 oil		●
9J-6(7)(8)	Stage plate (1) ; Stage plate (2); Stage plate (3)	●	●
9J-9	Stage plate (4)	○	○
9J.11	Stage plate holder	●	●
TBR03	Rotating binocular head (45 degree)	●	●
LH01	Lamp housing	●	●
	Halogen lamp 6V30W	●●	●●
PO02	Polarizer	●	●
PA01	Analyzer	●	●
CL02	Specimen clamps	●	●
9J.7	Photo adaptor	○	○
EPHp5	Photo eyepiece S5x	○	○
EPHp6.3	Photo eyepiece S6.3X	○	○
EPHp0.65	Video eyepiece 0.65X	○	○
	Metallurgical Microphotography equipments	○	○
	35mm DF camera with shutter release	○	○
SC03	Specimen holder	○	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139