

MIC-P5 Microscope

Specification		Model
		MIC-P5
Viewing Head	Compensation Free Trinocular Head Inclined at 30°(55mm-75mm)	●
Eyepiece	10×Gridding Eyepiece	●
	10×eyepiece with reticule	●
	WF10×/22mm	●
Nosepiece	Quintuple Nosepiece with centering adjustable	●
Objective	Stress-free Infinity Plan Objectives:4×,10×,25×,40×(s),60 ×(s)	●
Stage	Rotary Stage	●
	Stage Size: 172mm	●
	Moving Range: 360°	●
Condenser	Swing type N.A.1.25 Abbe condenser with iris diaphragm	●
Testing Piece	Gypsum 1λ testing piece	●
	Mica 1/4λ	●
	Quartz wedge testing piece	●
Focusing	Coaxial coarse & fine focusing adjustmen. With rack and pinion mechanism. Fine focusing scale value 0.002mm	●
Illuminator	High brightness fixed illumination	●
Light Source	Halogen Bulb 6V/30W.AC85V230V	●
	Adjustable Brightness	●
Checking Tool	0. 01mm Micrometer	●
Optional Accessory	Eyepiece: WF16×、WF20×、P16× 1.3Mega、2.0Mega、3.0Mega、 5.0Mega Pixels CMOS Digital camera eyepiece Stress-free objective:100×(s)Oil Photography attachment and CCD Adapter0.5×、×0.57×、0.75×	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

**CE ISO9001
ISO14001**



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139