

XT Serials Stereo Microscope

Model	XT-2A
Total Magnification	20X
Eyepieces	WF10X
Objectives	2X
Eyepieces Tube	Binocular Vertical Tube
Illumination	Natural Lamp
Coarse Adjustment Range	40mm
Working Distance	70mm
Working Stage	60mm Black/White, Plate and Clip



Model	XT-2B	XT-2C
Total Magnification	40X	40X
Eyepieces	WF10X	WF10X
Objectives	4X	4X
Eyepieces Tube	Binocular Vertical Tube	
Illumination	Natural Lamp	
Coarse Adjustment Range	40mm	40mm
Working Distance	60mm	60mm
Working Stage	60mm Black/White, Plate and Clip	

XT Serials Stereo Microscope

Model	XT-3B
Total Magnification	20X-40X
Eyepieces	WF10X
Objectives	2X 4X
Eyepieces Tube	Binocular 45° Inclined Tube with 55mm-75mm Distance Between Pupils
Illumination	Top light 12V/10W halogen lamp
Coarse Adjustment Range	40mm
Working Distance	57mm
Working Stage	60mm Black/White, Plate and Clip



Model	XT-3C
Total Magnification	20X-40X
Eyepieces	WF10X
Objectives	2X 4X
Eyepieces Tube	Binocular 45° Inclined Tube with 55mm-75mm Distance Between Pupils
Illumination	Top light 12V/10W halogen lamp
Coarse Adjustment Range	40mm
Working Distance	57mm
Working Stage	90mm frosted, glass or black/white, plate and clips

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

**CE ISO9001
ISO14001**



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139